



Accessories Catalogue

Disposable Range



www.neurogenmedical.com
info@neurogenmedical.com

About Us

Based in South Wales, Neurogen are the exclusive supplier of Deymed clinical neurophysiology equipment, to include Electroencephalography (EEG), Polysomnography (PSG), Electromyography (EMG), evoked potential (EP), ultrasound and transcranial magnetic stimulation (TMS) equipment. In addition, Neurogen offer a select range of neurophysiology consumables and accessories.

Neurogen have quickly built a reputation for having high-quality, industry-leading, products. This, paired with a focus on customer support, has led to high company growth with customer sites now spanning England, Wales, Scotland, Northern Ireland and the Republic of Ireland.

The Neurogen team have extensive experience working in clinical neurophysiology, neonatal intensive care, clinical research settings and the medical device industry. Company directors are qualified physiologists in Clinical Neurophysiology. They founded Neurogen as a platform to provide a direct service to neurophysiology colleagues across the UK and Ireland, prioritising customer service to support best in class patient care.

Our extensive range of disposable accessories are manufactured to the highest standards and are specifically designed to reduce potential cross-contamination.

**Neurogen**
www.neurogenmedical.com



Neurogen Limited
4 Castleton Court,
Fortran Road,
St Mellons, Cardiff,
South Wales
CF3 0LT

Tel: 029 2002 3316

Enquiries: info@neurogenmedical.com
Orders: orders@neurogenmedical.com

Table of Contents

EEG	5
Disposable EEG Cup Electrodes.....	6
Silver-Silver Chloride (Ag/AgCl).....	6
Standard Cup EEG Electrodes.....	6
Deep Cup EEG Electrodes.....	6
MRI and CT Compatible EEG Electrodes.....	7
Disposable Gold Cup Electrodes	8
Adult.....	8
Paediatric.....	8
Disposable EEG Application Kit.....	8
 EMG/EP	 9
Disposable NCS Electrodes	10
Single Pre-Gelled Adhesive Surface Electrodes	10
Twisted-Pair Pre-Gelled Adhesive Surface Electrodes	11
Disposable EMG Needle Electrodes	12
Technomed Concentric Needle Electrodes	12
MYOLINE Concentric Needle Electrodes	12
Disposable Macro and Single Fiber Needle Electrodes	13
Disposable Monopolar Needle Electrodes	13
Trocar Tip Needles.....	13
Pencil Tip Needles.....	14
Disposable Hypodermic Needle Electrodes.....	15
Disposable BAEP Accessories.....	15
Foam Ear Inserts	15
Disposable BAEP Headphone Covers.....	16
Disposable ERG Electrodes	16
ERG-JET Electrodes	16
ERG-Thread Electrodes	16
 PSG/SLEEP.....	 17
Nonin Elastic Sensor for SpO₂	18
Disposable Respiratory Effort Sensors	18
Disposable Airflow Thermal Sensor	18

IOM	19
Stimulation Probes.....	20
Monopolar Probes	20
Ball Tip Monopolar Probes.....	20
Bipolar Probes.....	21
Minifork	21
Ball Tip Bipolar Probes	22
Concentric Bipolar Probes.....	22
Hook Bipolar	23
Triple Hook Bipolar	23
DryTouch: Suction Stimulator Probe	23
Subdermal Needle Electrodes.....	24
Single Corkscrew.....	24
Twisted Corkscrew.....	24
Disposable Subdermal Hooked.....	25
Disposable Subdermal Curved.....	25
Subdermal Needle Electrodes	25
Single Subdermal Needle Electrodes.....	25
MRI Single Subdermal Needle Electrodes	26
Paired Subdermal Needles.....	26
3-Wire Twisted Subdermal Needles.....	26
5-Wire Twisted Subdermal Needles.....	27
Laryngeal Electrodes.....	27

EEG

TruScan® EEG
EEG / EP / LTM | SYSTEM



- Compact headbox with wireless and ambulatory options
- Dedicated EEG Review Keyboard
- Optical Isolation for industry-leading signal quality
- Intuitive "Dashboard" interface
- Always-On Impedance Indication
- Up to 3KHz sampling rate
- Up to 256-channels
- Routine, LTM, Portable, Wireless, ICU, Ambulatory and Video Ambulatory/HVT configurations

Disposable EEG Cup Electrodes

User Benefits



- No cleaning and preparation of EEG electrodes needed
- Clear and reliable signal due to optimised layer of frequency stable Ag/AgCl
- Optimised signal quality due to larger electrode surface
- Short cup-to-cable connection to reduce lift risk
- Optimised connector for easy headbox insertion and extraction
- Optimal patient comfort with soft and flexible anti-tangle lead wires
- Reduced risk of cross-contamination
- Latex-free, Lead-free and BPA-free
- REACH and RoHS Compliant
- Standard DIN 42802 Touch-Proof Connectors

Silver-Silver Chloride (Ag/AgCl)

Standard Cup EEG Electrodes

Part Number	Lead Length (cm)	Cup Diameter (mm)	QTY
TE/C52-432	100	10	10 pcs/pack
			100 pcs/box
TE/C52-632	150	10	10 pcs/pack
			100 pcs/box
TE/C52-832	200	10	10 pcs/pack
			100 pcs/box
TE/C52-932	250	10	10 pcs/pack
			100 pcs/box
TE/C52-438	100	10	25 pcs/pack
			250 pcs/box
TE/C52-638	150	10	25 pcs/pack
			250 pcs/box
TE/C52-838	200	10	25 pcs/pack
			250 pcs/box
TE/C52-938	250	10	25 pcs/pack
			250 pcs/box



Part Number	Lead Length (cm)	Cup Diameter (mm)	QTY
DAGD102600	100	10	10 pcs
DAGD152600	150	10	10 pcs
DAGD252600	250	10	10 pcs
DAGD102600/25	100	10	25 pcs
DAGD152600/25	150	10	25 pcs
UNIMED-DISP-EEG10	100	10	25 pcs
UNIMED-DISP-EEG15	150	10	25 pcs



Deep Cup EEG Electrodes

- Specifically designed to hold a higher volume of conductive paste.



Part Number	Lead Length (cm)	Cup Diameter (mm)	QTY
TE/C53-432	100	10	10 pcs/pack
			100 pcs/box
TE/C53-632	150	10	10 pcs/pack
			100 pcs/box
TE/C53-832	200	10	10 pcs/pack
			100 pcs/box
TE/C53-932	250	10	10 pcs/pack
			100 pcs/box
TE/C53-438	100	10	25 pcs/pack
			250 pcs/box
TE/C53-638	150	10	25 pcs/pack
			250 pcs/box
TE/C53-838	200	10	25 pcs/pack
			250 pcs/box
TE/C53-938	250	10	25 pcs/pack
			250 pcs/box

MRI and CT Compatible EEG Electrodes

- Specifically designed to use in MRI or CT diagnostic techniques.

Part Number	Technique	Lead Length (cm)	Cup Diameter (mm)	QTY
DAGR102600MRI	MRI	100	10	10 pcs
DAGR152600MRI	MRI	150	10	10 pcs
DAGR202600MRI	MRI	200	10	10 pcs
DAGR102600	CT	100	10	10 pcs
DAGR152600	CT	150	10	10 pcs
DAGR202600	CT	200	10	10 pcs



MRI, DAGRXX2600MRI



Disposable Gold Cup Electrodes



Adult

- Moulded cable-electrode connection ensures high mechanical strength during use.

Part Number	Lead Length (cm)	Cup Diameter (mm)	QTY
DEGM052600	50	10	10 pcs
DEGM102600	100	10	10 pcs
DEGM152600	150	10	10 pcs
DEGM252600	250	10	10 pcs

Paediatric

- Designed for paediatric patients.
- Reduced diameter of electrode to allow recording in more focused areas.

Part Number	Lead Length (cm)	Cup Diameter (mm)	QTY
DEGS052606	50	6	10 pcs
DEGS102606	100	6	10 pcs
DEGS152606	150	6	10 pcs
DEGS252606	250	6	10 pcs


Disposable EEG Application Kit

- EEG Application kits, designed to significantly reduce skin breakdown, cross contamination, and excess waste.

Part Number	Description	QTY
KITECH001	TechSystem Kit, includes a set of 24 Tech Dots, a single use Tech Prep stick, and 24 individual Tech Tape pads for added adhesive	1
DOT024-25	Tech Dots (24 dots/set)	25 pcs
NPREP003K-25	Tech Prep stick	25 pcs



KITECH001



DOT024-25

EMG/EP

TruTrace® EMG NCS / EMG / EP | SYSTEM



- Modern design with intelligent charging battery operation for superior signal quality
- Fast and easy to use with dashboard style interface
- Basic to advanced configuration of all tests
- Full control of all functions via dedicated EMG Keyboard

NCS

EMG

Quantitative EMG

Multi EMG

Macro EMG

SF EMG

F-Wave

Ultrasound

SSEP

MEP

VEP

BAEP

Cognitive

MUNE/MUNIX

Autonomic

Desktop EMG/EP & Ultrasound



Disposable NCS Electrodes

Single Pre-Gelled Adhesive Surface Electrodes

1.5mm Touch-Proof Connectors (42802)

Part Number	Size (mm)	Lead length (cm)	QTY
Recording Electrodes			
TE-K50430-001	20	100	4 pcs/pack 40 pcs/box
TE-K50430-015	20	150	4 pcs/pack 40 pcs/box
TE/K50430-002	20	200	4 pcs/pack 40 pcs/box
TE-K50430-025	20	250	4 pcs/pack 40 pcs/box
DENIS10026	15 x 20	100	20 pcs/pack 120 pcs/box
DENIS15026	15 x 20	150	20 pcs/pack 120 pcs/box
DENIS25026	15 x 20	250	20 pcs/pack 120 pcs/box
DENIL10026	20 x 25	100	20 pcs/pack 120 pcs/box
DENIL15026	20 x 25	150	20 pcs/pack 120 pcs/box
DENIL25026	20 x 25	250	20 pcs/pack 120 pcs/box
Ground Electrodes			
TE/K50431-002	40 x 50	100	1 pc/pack 20 pcs/box
TE/K50431-003	40 x 50	200	1 pc/pack 20 pcs/box



0.7mm Touch-Proof Connectors

Part Number	Description	Lead length (cm)	QTY
97-153	Pre-gelled Ag-AgCl, Friendship Medical	10	12pcs/pack 120 pcs/box
97-000	Pre-gelled Ag-AgCl, Dantec	10	12pcs/pack 120 pcs/box


Reusable connection cable (97-105 or 97-100) required



Tab and Snap Adhesive Surface Electrodes

- For use with Crocodile Clip and Snap connectors.

Part Number	Size (mm)	Connectors	QTY
TABSMPO00025	14 X 34	Tab	100 pcs/pack 2500 pcs/box
TABSM000025	22 X 34	Tab	100 pcs/pack 2500 pcs/box
CDES000024	24	Snap	25 pcs/pack 500 pcs/box



CDES000024

Twisted-Pair Pre-Gelled Adhesive Surface Electrodes

Standard Touch-Proof connectors

Part Number	Size (mm)	Lead length (cm)	QTY
96-238	15 x 20	2x 150	24 pcs
96-239	15 x 20	2x 250	24 pcs
96-240	20 x 27	2x 150	24 pcs
96-241	20 x 27	2x 250	24 pcs

Disposable EMG Needle Electrodes

Technomed Concentric Needle Electrodes

Part Number	Needle Length (mm)	Diameter (mm)	Rec area (mm ²)	Colour	QTY
TE/B50600-002	25	0.40 (27G)	0.031		25 pcs
TE/B50600-001	25	0.30 (30G)	0.031		25 pcs
TE/B50600-006	30	0.30 (30G)	0.031		25 pcs
TE/B50600-003	37	0.45 (26G)	0.085		25 pcs
TE/B50600-004	50	0.45 (26G)	0.085		25 pcs
TE/B50600-005	75	0.60 (23G)	0.085		25 pcs

Reusable connection needle cable (TE/D64-551 or TE/D64-851) required

P/N: TE/D64-551, 1.2m

P/N: TE/D64-851, 2m



MYOLINE Concentric Needle Electrodes

Part Number	Needle Length (mm)	Diameter (mm)	Rec area (mm ²)	Colour	QTY
D0390302511	25	0.30 (30G)	0.02		25 pcs
D039035302	30	0.35 (28G)	0.07		25 pcs
D0390453010	30	0.45 (26G)	0.07		25 pcs
D039035408	40	0.35 (28G)	0.07		25 pcs
D039045401	40	0.45 (26G)	0.07		25 pcs
D039045504	50	0.45 (26G)	0.07		25 pcs
D039060706	75	0.60 (23G)	0.07		25 pcs

Reusable connection needle cable (100X00NNHF3 or 100200NNHF3) required


P/N: 10000NNHF3, 1m

P/N: 100200NNHF3, 2m



Disposable Macro and Single Fiber Needle Electrodes

- Recording area on a small side opposite to the cutting edge and close from the tip.
- The Macro version is made with insulated cannula.



Part Number	Needle Length (mm)	Diameter (mm)	Rec area (mm²)	Model	QTY
DSFN380452506	25	0.45 (26G)	0.0011	Single Fiber	5 pcs
DSFN380453707	37	0.45 (26G)	0.0011	Single Fiber	5 pcs
Reusable connection needle cable (10000NNHF3 or 100200NNHF3) required for use with SF needle					
DSFN380453705	37	0.45 (26G)	0.0011	Macro	5 pcs
Reusable connection cable (100000NNHFM, 1m) required for use with Macro needle					

P/N: 10000NNHF3, 1m
P/N: 100200NNHF3, 2m



P/N: 10000NNHFM, 1m



Disposable Monopolar Needle Electrodes

Trocars Tip Needles



Pre-wired

Cable length 75cm



Part Number	Needle Length (mm)	Diameter (mm)	Colour	QTY
TE/AP2535-335/01	25	0.35 (28G)	●	25 pcs
TE/AP2545-335/01	25	0.45 (26G)	●	25 pcs
TE/AP3735-335/01	37	0.35 (28G)	●	25 pcs
TE/AP3745-335/01	37	0.45 (26G)	●	25 pcs
TE/AP5045-335/01	50	0.45 (26G)	●	25 pcs
TE/AP7545-335/01	75	0.45 (26G)	●	25 pcs

Detachable



Part Number	Needle Length (mm)	Diameter (mm)	Colour	QTY
TE/AD2535-008/01	25	0.35 (28G)	●	25 pcs
TE/AD2545-008/01	25	0.45 (26G)	●	25 pcs
TE/AD3735-008/01	37	0.35 (28G)	●	25 pcs
TE/AD3745-008/01	37	0.45 (26G)	●	25 pcs
TE/AD5045-008/01	50	0.45 (26G)	●	25 pcs
TE/AD7545-008/01	75	0.45 (26G)	●	25 pcs

Detachable needles require monopolar needle cable, P/N: TE/D54-451 or TE/D34-231




TE/D54-451, 1m



TE/D34-231, 0.6m

Pencil Tip Needles



Pre-wired

Cable length 75cm



Part Number	Needle Length (mm)	Diameter (mm)	Colour	QTY
TE/AP2535-335	25	0.35(28G)	●	25 pcs
TE/AP2545-335	25	0.45 (26G)	●	25 pcs
TE/AP3735-335	37	0.35 (28G)	●	25 pcs
TE/AP3745-335	37	0.45 (26G)	●	25 pcs
TE/AP5045-335	50	0.45 (26G)	●	25 pcs
TE/AP7545-335	75	0.45 (26G)	●	25 pcs

Detachable



Part Number	Needle Length (mm)	Diameter (mm)	Colour	QTY
TE/AD2535-008	25	0.35 (28G)	●	25 pcs
TE/AD2545-008	25	0.45 (26G)	●	25 pcs
TE/AD3735-008	37	0.35 (28G)	●	25 pcs
TE/AD3745-008	37	0.45 (26G)	●	25 pcs
TE/AD5045-008	50	0.45 (26G)	●	25 pcs
TE/AD7545-008	75	0.45 (26G)	●	25 pcs

Detachable needles require monopolar needle cable, P/N: TE/D54-451 or TE/D34-231



TE/D54-451, 1m



TE/D34-231, 0.6m

Disposable Hypodermic Needle Electrodes

For application of Botulinum Toxin

User Benefits

- Optimal patient comfort with ultra-sharp tip and special coating for very smooth penetration and repositioning.
- Leak-free and minimum drug residual volume.
- Optimised to fit to all slip-tip and Luer Lock syringe connections.

Part Number	Needle Length (mm)	Diameter (mm)	Colour	QTY
TE/L2530-335	25	0.30 (30G)		10 pcs
TE/L2540-335	25	0.40 (27G)		10 pcs
TE/L3740-335	37	0.40 (27G)		10 pcs
TE/L3745-335	37	0.45 (26G)		10 pcs
TE/L5050-335	50	0.50 (25G)		10 pcs
TE/L7570-335	75	0.70 (22G)		10 pcs

Part Number	Needle Length (mm)	Diameter (mm)	Colour	QTY
MIOBOT2530	25	0.30 (30G)		10 pcs
MIOBOT2540	25	0.40 (27G)		10 pcs
MIOBOT3745	37	0.40 (27G)		10 pcs
MIOBOT5050	50	0.45 (26G)		10 pcs
MIOBOT7570	75	0.70 (22G)		10 pcs

Disposable BAEP Accessories

Foam Ear Inserts

Part Number	Description	QTY
1700-9636-500D	E-A-RLINK™ 3B Foam Insert Eartips, Beige Small Tips	500 pcs
1700-9636-750D	E-A-RLINK™ 3B Foam Insert Eartips, Beige Small Tips	750 pcs
1700-9635-500D	E-A-RLINK™ 3A Foam Insert Eartips, Yellow Standard Tips	500 pcs
1700-9635-750D	E-A-RLINK™ 3A Foam Insert Eartips, Yellow Standard Tips	750 pcs



Disposable BAEP Headphone Covers

Part Number	Description	QTY
HYGENX50	HygenX Sanitary Ear Cushion Covers for Over-Ear Headphone	50 pcs



Disposable ERG Electrodes ERG-JET Electrodes

Part Number	Description	Lead Length (cm)	QTY
F-06 CHF-12	Gold foil electrode on poly (methyl methacrylate) with Touch-Proof connector, 7.9 Radius of curve	100	12 pcs
F-06 CHF-50	Gold foil electrode on poly (methyl methacrylate) with Touch-Proof connector, 7.9 Radius of curve	100	50 pcs



F-06 CHF-12


ERG-Thread Electrodes

Part Number	Description	Lead Length (cm)	QTY
DT0001	GVB-Spes Single use sterile ERG corneal thread electrode, GVB	15	10
ET-01	Fabrinal Single use sterile ERG corneal thread electrode. <i>This does not require an extension cable</i>	100	10

DT0001 requires TP extension cable



ET-01



DT0001

SomniPro® PSG

PSG / EEG / LTM | SYSTEM

- Small-lightweight Wireless PSG and Full EEG
- Wireless mode option with complete synchronisation to Full-HD Video
- Long-lasting battery operation, even in wireless mode
- Height-adjustable cart
- Passive cooling, silent, PC



PSG/Sleep

Disposable Sleep Accessories

Nonin Elastic Sensor for SpO₂

- Disposable PureLight sensors.

Part Number	Description	QTY
95-600	Adult Sensor	24 pcs
95-610	Paediatric Sensor	24 pcs
95-755	Adult Flexi-wrap	25 pcs
95-756	Paediatric Flexi-wrap	25 pcs



Disposable Respiratory Effort Sensors

Part Number	Description	Lead Length (cm)	QTY
RESPA000200	Adult	200	5 pcs
RESPP000200	Paediatric	200	5 pcs

Requires Connection Cable BMDS00260R, 60cm cable length



RESPP000200



RESPA000200




BMDS00260R


Disposable Airflow Thermal Sensor

- Breathing sensor based on the thermal principle.

Part Number	Description	Lead Length (cm)	QTY
AIRFA00300	Airflow - Adult	200	1 pc
AIRFAP0300	Airflow - Infant	200	1 pc



AIRFA00300



AIRFAP0300

IOM Accessories



IOM

Stimulation Probes


Monopolar Probes

MP100/1



Part Number	Model	Probe Length (mm)	Diameter (mm)	QTY
MP100/1	Monopolar	100	0.75	1 pc
MP100/10	Monopolar	100	0.75	10 pcs
MP200/1	Monopolar	200	1.5	1 pc
MP200/10	Monopolar	200	1.5	10 pcs

Part Number	Model	Probe Length (mm)	QTY
SI1I0045S2526D	Monopolar	45	5 pcs
SI1I0100S2526D	Monopolar	100	5 pcs
SI1I0200S2526D	Monopolar	200	5 pcs
SI1IB130S2526D	Bayonet monopolar	130	5 pcs
SI1IC095S2526D	Bent monopolar	95	5 pcs




Ball Tip Monopolar Probes

BTP2.3-130/1



Part Number	Model	Probe Length (mm)	Diameter (mm)	QTY
BTP3.0-100/1	Monopolar	100	3	1 pc
BTP3.0-100/10	Monopolar	100	3	10 pcs
BTP3.0-200/1	Monopolar	200	3	1 pc
BTP3.0-200/10	Monopolar	200	3	10 pcs
BTP2.3-130/1	Monopolar	130	2.3	1 pc
BTP2.3-130/10	Monopolar	130	2.3	10 pcs
BTP2.3-200/1	Monopolar	200	2.3	1 pc
BTP2.3-200/10	Monopolar	200	2.3	10 pcs

Part Number	Model	Probe Length (mm)	Diameter (mm)	QTY
SI1B0045S2526D	Ball Monopolar	45	2	5 pcs
SI1B0100S2526D	Ball Monopolar	100	3	5 pcs
SI1BB130S2526D	Ball Bayonet Monopolar	130	3	5 pcs



Bipolar Probes

Minifork



BIP100/1

Part Number	Model	Shaft Length (mm)	Diameter (mm)	QTY
BIP100/1	Bipolar Minifork	100	0.75	1 pc
BIP100/10	Bipolar Minifork	100	0.75	10 pcs

Part Number	Model	Shaft Length (mm)	QTY
SI2F0045S2526D	Bipolar Minifork	45	5 pcs
SI2F0100S2526D	Bipolar Minifork	100	5 pcs
SI2FB130S2526D	Bayonet Bipolar Minifork	130	5 pcs
SI2FE095S2526D	90° Bent Bipolar Minifork	95	5 pcs
SI3FE095S2526D	90° Bent Tripolar Paediatric Minifork	95	5 pcs




Ball Tip Bipolar Probes

Part Number	Model	Probe Length (mm)	Diameter (mm)	QTY
SI2B0100S2526D	Ball Bipolar	100	3	5 pcs
SI2B0030S2526D	Ball Bipolar	30	2	5 pcs




Concentric Bipolar Probes



Part Number	Model	Probe Length (mm)	Diameter (mm)	QTY
CP100/1	Concentric Bipolar	100	1	1 pc
CP100/10	Concentric Bipolar	100	1	10 pcs

Part Number	Model	Probe Length (mm)	QTY
SI2C0045S2526D	Concentric Bipolar	45	5 pcs
SI2C0100S2526D	Concentric Bipolar	100	5 pcs
SI2CB130S2526D	Bayonet Concentric Bipolar	130	5 pcs
SI2CC045S2526D	Bent Concentric Bipolar	45	5 pcs
SI2CC095S2526D	Bent Concentric Bipolar	95	5 pcs
SI2TC045S2526D	Ball Concentric Bipolar	45	5 pcs
SI2TC100S2526D	Ball Concentric Bipolar	100	5 pcs



Part Number	Model	Shaft Length (mm)	QTY
SI1IF095S2526D	Flush Tip	95	5 pcs



Hook Bipolar

Part Number	Model	Shaft Length (mm)	QTY
SI2H0020S2526D	Hook	20	5 pcs



Triple Hook Bipolar

Part Number	Model	Shaft Length (mm)	QTY
SI3H0020S2526D	Triple Hook	20	5 pcs



DryTouch: Suction Stimulator Probe

- Provides a dry stimulation field for accurate testing.
- Tip designed to simultaneously provide repeatable stimulation and evacuate conductive fluids that generate electrical shunting.

Part Number	Model	Shaft Length (mm)	QTY
SI1IK0260S2526D	Dry Touch	260	5 pcs
SI1IK0130S2526D	Dry Touch	130	5 pcs



Subdermal Needle Electrodes

Single Corkscrew

Part Number	Description	Lead Length (cm)	Colour	QTY
TE/S50715-001	Singles	120	6 colour combinations (Red, Green, Black, White, Yellow, Blue)	24 pcs
MNCS0615DS	Singles	150	6 colour combinations (Red, Green, Black, White, Yellow, Blue)	24 pcs
TE/S50715-002	Singles	200	6 colour combinations (Red, Green, Black, White, Yellow, Blue)	24 pcs
TE/S50715-003	Singles	250	6 colour combinations (Red, Green, Black, White, Yellow, Blue)	24 pcs


MRI (Magnetic Resonance Imaging) compatible version available on request



Twisted Corkscrew

Part Number	Description	Lead Length (cm)	Colour	QTY
96-237	Twisted x2	2x 200	12 colour combinations	12 pcs
MNCS0610DTS	Twisted x2	2x 100	12 colour combinations	12 pcs
MNCS0615DTS	Twisted x2	2x 150	12 colour combinations	12 pcs
MNCS0625DTS	Twisted x2	2x 250	12 colour combinations	12 pcs
MNCS0615DT5	Twisted x5	5x 150	5 colour combinations	10 pcs
MNCS0615DT6	Twisted x6	6x 150	6 colour combinations	10 pcs


MRI (Magnetic Resonance Imaging) compatible version available on request



Disposable Subdermal Hooked

4 x 6 colours (Red, Green, Black, White, Yellow, Blue)


Part Number	Description	Lead Length (cm)	QTY
TE/S50719-001	Singles	150	24 pcs
TE/S50719-002	Singles	200	24 pcs
TE/S50715-003	Singles	250	24 pcs



Disposable Subdermal Curved

5 x 5 colours (Red, Green, White, Yellow, Blue)


Part Number	Description	Lead Length (cm)	QTY
TE/S45-638	Singles	150	25 pcs
TE/S45-838	Singles	200	25 pcs



Subdermal Needle Electrodes

Single Subdermal Needle Electrodes


4 x 6 colours (Red, Green, Black, White, Yellow, Blue)



Part Number	Needle Length (mm)	Diameter (mm)	Lead Length (cm)	QTY
TE/S50716-001	14	0.38	100	24 pcs
TE/S50716-002	14	0.38	150	24 pcs
TE/S50716-003	14	0.38	200	24 pcs
TE/S50716-004	14	0.38	250	24 pcs
TE/S43-438	13	0.4	100	24 pcs
TE/S43-638	13	0.4	150	24 pcs
TE/S43-838	13	0.4	200	24 pcs
TE/S43-938	13	0.4	250	24 pcs

MRI Single Subdermal Needle Electrodes

Part Number	Needle Length (mm)	Diameter (mm)	Lead Length (cm)	QTY
MN4013D10SMRI	13	0.4	100	24 pcs
MN4013D15SMRI	13	0.4	150	24 pcs
MN4013D25SMRI	13	0.4	250	24 pcs
MN4018D10SMRI	18	0.4	100	24 pcs
MN4018D15SMRI	18	0.4	150	24 pcs
MN4018D25SMRI	18	0.4	250	24 pcs



Paired Subdermal Needles



Part Number	Description	Needle Length (mm)	Diameter (mm)	Lead Length (cm)	QTY
TE/S44-637	Twisted	13	0.4	2x 150	20 pcs
TE/S44-837	Twisted	13	0.4	2x 200	20 pcs
TE/S44-937	Twisted	13	0.4	2x 250	20 pcs
TE/S22-837	Twisted	22	0.3	2x 200	8 pcs
TE/S48-637	Parallel	13	0.4	2x 150	10 pcs
TE/S48-837	Parallel	13	0.4	2x 200	10 pcs
TE/S48-937	Parallel	13	0.4	2x 250	10 pcs

Available as a colour mix


3-Wire Twisted Subdermal Needles



Part Number	Needle Length (mm)	Diameter (mm)	Lead Length (cm)	QTY
TT03MN4013D/10	13	0.4	3x 100	10 pcs
TT03MN4013D/15	13	0.4	3x 150	10 pcs
TT03MN4013D/25	13	0.4	3x 250	10 pcs
TT03MN3518D/10	18	0.35	3x 100	10 pcs
TT03MN3518D/15	18	0.35	3x 150	10 pcs
TT03MN3518D/25	18	0.35	3x 250	10 pcs

Available as a colour mix

5-Wire Twisted Subdermal Needles



Part Number	Needle Length (mm)	Diameter (mm)	Lead Length (cm)	QTY
TT05MN4013D/10	13	0.4	5x 100	10 pcs
TT05MN4013D/15	13	0.4	5x 150	10 pcs
TT05MN4013D/25	13	0.4	5x 250	10 pcs
TT05MN3518D/10	18	0.35	5x 100	10 pcs
TT05MN3518D/15	18	0.35	5x 150	10 pcs
TT05MN3518D/25	18	0.35	5x 250	10 pcs

Available as a colour mix

Laryngeal Electrodes

Dragonfly – Adhesive Laryngeal Electrodes

- The Dragonfly electrode is made of soft material to avoid any risk of injury while inserting the tube.
- Can be applied on any reinforced endotracheal tube.
- 1.5mm Touch Proof Connectors (42802).

Part Number	Number of channels	ET Tube Diameter (mm)	Lead Length (cm)	Pack	
LSE3446MCP075	1	2.0 – 5.5	75	10 pcs	
LSE3446MCA075	1	6.0 – 7.0	75	10 pcs	
LSE4246MCA075	1	7.5 – 10.0	75	10 pcs	
LSE3463DCA200	2	6.0 – 7.5	200	10 pcs	
LSE4263DCA200	2	8.0 – 9.5	200	10 pcs	


Integrated Laryngeal Electrodes

- 1.5mm Touch Proof Connectors (42802).
- 200cm cable length.



Armed Version

Part Number	Number of channels	ET Tube Diameter (mm)	Pack
LSEPA01	1	6	5 pcs
LSEPA02	1	7	5 pcs
LSEPA03	1	8	5 pcs
LSEPA05	2	6	5 pcs
LSEPA06	2	7	5 pcs
LSEPA07	2	8	5 pcs



Not Armed Version

Part Number	Number of channels	ET Tube Diameter (mm)	Pack
LSEP001	1	6	5 pcs
LSEP002	1	7	5 pcs
LSEP003	1	8	5 pcs
LSEP005	2	6	5 pcs
LSEP006	2	7	5 pcs
LSEP007	2	8	5 pcs



Neurogen Limited
4 Castleton Court,
Fortran Road, St Mellons
Cardiff, South Wales
CF3 0LT

Tel: 029 2002 3316

Email: info@neurogenmedical.com